

# Genepi Molmed Hackathon

---

Welcome to the Genepi MinION Hackathon! Your task today is to analyse mtDNA data generated from Oxford Nanopore. The data has been generated using a transposase library prep kit on a mixture of two long-range PCR amplicons, which each amplify about half of the mtDNA genome. The amplicons overlap by about 200 at each end. Mind that a transposase library prep on so short amplicons (about 8.5 kb each) produces rather short reads. These are rather untypical for nanopore sequencing at genome scale or when using a ligation library prep kit, which would specifically ligate the adaptors to the end of the amplicons and thus produce full-length amplicon reads.

## Setup

- Create a Word Document with your last name (named: `SS24_NGS_Molmed_Lastname.docx`) and document what you do (e.g. commands, screenshots, graphics). This file must be sent at the end of the class to Sebastian.
- Login to `genepi-lehre.i-med.ac.at` and create a folder `project-day` within `genepi-teaching/students/<your-q-number>`. All experiments must be executed from here.
- Copy your barcode folder (`cp -r <your-barcode> .`) to your `project-day` folder. The files are located here: `~/genepi-teaching/ngs/2024-ngs-molmed`. The structure should look like below (e.g. for barcode04).
- After the `cp` has been finished, verify that fastq files are available within your folder.

## Structure

```
└─ project-day
  └─ barcode04
    └─ file1.fastq.gz
    └─ file2.fastq.gz
    └─ ...
```

## Quality Control

### Task 1

Count the reads in your sequencing experiment by combining the commands `zcat *fastq.gz |`, and `wc -l` and explain the command and output in your Word file. If you see almost no reads included (< 400), copy a different barcode to your `project-day` folder and redo Task 1.

### Task 2

The second task is to run a QC analysis on your data. For Nanopore data, we use a tool called [NanoPlot](#). The usage looks complicated but at the end, you only need the following options: `--fastq-rich <fastq-files>, --outdir <output-folder>, --maxlength 15000`).

```
NanoPlot --fastq-rich <your-fastq-files> ...
```

**HINT:** Use `*fastq.gz` if you have several files!

This tool outputs a HTML report, look for the report, download it to your local filesystem and see what information you can identify. Add some plots and information to your Word file!

## Mapping

Now it's time to align the data using a tool called [minimap2](#). Replace with your input-file (included in your barcode folder) and with the name of your actual barcode. The reference can be found here: [~/genepi-teaching/students/reference-data/chrM.fasta](#)

```
minimap2 -ax map-ont <reference> <your-fastq-files> | samtools sort >
<yourbarcode>.bam
```

**HINT:** Use `*fastq.gz` for `<your-fastq-files>` if you have several files!

Copy the command into your doc-file and describe what the command is doing (e.g. `minimap2`, `|`, `samtools sort`).

## Visualize Mapping

Below we execute a new command that calculates the coverage for each position. Run `samtools depth` on your BAM file. The output is written on the command-line by default. Adapt the command and write the result into a file (using `>`), transfer the file to Windows and visualize the result.txt file with Excel or R.

```
samtools depth -a <yourbamfile>
```

## Variant Calling

Now we run the actual variant calling with a tool created specific for mitochondrial genomes.

```
/opt/tools/mutserve/mutserve call <your-bam.bam> --reference <reference> -
-level 0.1 --output <output-name>.vcf.gz
```

Questions:

- What is `--level 0.1` doing? Try different parameters.
- Open the vcf file with `zcat` and see what information you can identify. Discuss it in your Word file. The VCF file format is explained in the slides or also [here](#).

## Haplogroups

mtDNA profiles can also be grouped into so called [haplogroups](#). Download the file to your local PC and upload the VCF File to <https://haplogrep.i-med.ac.at>.

Questions:

- Interpret the results (use the Haplogrep docs) and add information to your Word file.

PCR Products

